IVF Basics® SpermTec M

IVF Basics® SpermTec M is a ready-to-use HEPES-buffered medium which contains bicarbonate, physiologic salts, glucose, lactate, pyruvate, phenol red, 10mg/liter gentamicin sulphate (medicinal substance) and 5g/liter HSA (medicinal substance derived from human blood plasma).

The MSDS for IVF Basics® SpermTec M also describes the product composition and is available on the website of Gynotec B.V.

IVF Basics® SpermTec M has been developed for all human sperm washing and sperm preparation techniques.

IVF Basics® SpermTec M contains HEPES, no CO₂ incubation is required.

APPLICATION

IVF Basics® SpermTec M is designed for all human sperm washing and sperm preparation techniques. It can also be used for intrauterine insemination (IUI).

For professional use only.

COMPOSITION

IVF Basics® SpermTec M is a ready-to-use HEPES-buffered medium which contains bicarbonate, physiologic salts, glucose, lactate, pyruvate, phenol red, 10mg/liter gentamicin sulphate (medicinal substance) and 5g/liter HSA (medicinal substance derived from human blood plasma).

The MSDS for IVF Basics® SpermTec M also describes the product composition and is available on the website of Gynotec B.V.

QUALITY CONTROL

- pH: 7.20-7.60
- Osmolality: 270-290 mOsm/kg
- Endotoxins test by LAL methodology (USP <85>): < 0.25 EU/ml
- Sterility test by the current Ph. Eur. 2.6.1/ USP <71>: No growth
- One-cell mouse embryo assay (% blastocysts after 96h): ≥ 80%
- Human sperm survival assay (% motility compared with control after 24 hours): ≥ 80%
- Chemical composition
- Use of Ph Eur or USP grade products if applicable
- Certificate of analysis and MSDS are available upon request

Sterilized by using aseptic processing techniques.

PRECAUTIONS AND WARNINGS

- Aseptic techniques should be used to avoid possible contamination, even when the products contain gentamicin. Always wear protective clothing when handling specimens.
- Handle specimens as if capable of transmitting HIV or hepatitis.
- All blood products should be treated as potentially infectious. Source material used to manufacture this product was tested and found non-reactive for HbsAg and negative for Anti-HIV-1/-2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents.
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of virus. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specification by established processes.
- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to Gynotec B.V. and, if applicable, to the competent authority of the EU Member State in which the user and/or patient is established.

PRE-USE CHECKS

- Do not use the product if:
 - the seal on the bottle is broken or open when the product is delivered;
 - it becomes discoloured, shows any evidence of microbial contamination or cloudy;
 - expiry date has been exceeded.
- Do not freeze before use.
- Do not re-sterilize after opening.
- Products which include gentamicin should not be used on a patient that has a known allergy to gentamicin or similar antibiotics.
- Depending on the number of procedures that will be performed on one day, remove the required volume of medium under aseptic conditions in an appropriate sterile recipient. This is in order to avoid multiple openings/warming cycles of the medium. Discard excess (unused) media.

STORAGE CONDITIONS

- Store refrigerated (2°C – 8°C).
- Keep away from (sun)light.
- After opening the container, do not use the product longer than 7 days. Sterile conditions must be maintained and product must be stored at 2°C – 8°C.
- Discard the devices in accordance with local regulations for disposal of medical devices.

IVF Basics® SpermTec M

ID: FP09 I70 RO1_Sperm C.1
Update: 01/08/2023

IVF Basics® SpermTec M contains a HEPES buffer to maintain stable pH during procedures performed outside a CO₂ incubator. Incubation in a CO₂ incubator will lower the pH.

INSTRUCTIONS FOR USE

Each laboratory should establish its own validated procedures and protocols.

IVF Basics® SpermTec M contains a HEPES buffer to maintain stable pH during procedures performed outside a CO₂ incubator. Incubation in a CO₂ incubator will lower the pH.

Density gradient centrifugation

- Bring all components to room temperature or 37°C.
- Centrifuge at 2100 rpm (700g) for 15 minutes.
- Carefully aspirate the supernatant from the pellet.
- Open the IVF Basics® SpermTec M.
- Add with a pipette 2-4ml IVF Basics® SpermTec M on top and resuspend the pellet by shaking.
- Centrifuge at 1500 rpm (400g) for 7 minutes.
- Centrifuge at 2100 rpm (700g) for the pellet, leaving the pellet intact.
- Add 6-4ml IVF Basics® SpermTec M and resuspend. Perform a microscopic analysis.
- Store the tube at room temperature or in an incubator protected against sun light until just before insemination.
- Before insemination centrifuge at 1100 rpm (200g) for 7 minutes.
- Carefully aspirate the supernatant until a volume of 0.2-0.5ml above the pellet.
- Resuspend the pellet by shaking.
- Ready for insemination.

Swim-up (According to WHO, 2021)

- Carefully layer 1.5ml IVF Basics® SpermTec M on top of 1ml washed semen in a conical bacteriological centrifuge tube.
- Place the tube in the incubator at a 45° angle for 1 hour at 37°C.
- Place the tube carefully in an upright position and remove 1ml of the top layer.
- Dilute the aliquot of motile cells with 2-5ml of IVF Basics® SpermTec M. Centrifuge for 5-10 minutes at 300g.
- Carefully aspirate the supernatant and resuspend the sperm pellet in 0,5ml IVF Basics® SpermTec M.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)
The SSCP for IVF Basics® SpermTec M describes safety and performance characteristics for the medium, and is available on the website of Gynotec B.V. (www.gynotec.nl).

For further questions regarding to the safety and performance, please contact Gynotec B.V.

IVF Basics® SpermTec M

ID: FP09 I70 RO1_Sperm C.1
Update: 01/08/2023

IVF Basics® SpermTec M contains a HEPES buffer to maintain stable pH during procedures performed outside a CO₂ incubator. Incubation in a CO₂ incubator will lower the pH.

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IVF Basics® SpermTec M için MSDS aynı zamanda ürün bileşimini de açıklar ve Gynotec B.V.'nin web sitesinde mevcuttur.