	SpermWash®
FP09 I08_SW R01 A.1	
Date: 25/07/2023	
APPLICATION	
SpermWash® is intended for in vitro procedures involving human gametes (sperm and oocytes), including washing of gametes, sperm swim-up procedures, intra-uterine insemination (IUI) of the spermatozoa and sperm injection during intracytoplasmic sperm injection (ICSI). SpermWash® can also be used for human embryo washing and holding, and for embryo transfer (ET) in the uterine cavity.	
<i>For professional use only.</i>	
COMPOSITION	
SpermWash® is a ready-to-use HEPES buffered medium which also contains phenol red, bicarbonate, physiologic salts, glucose, lactate, pyruvate and Human Serum Albumin (4.0 g/L; medicinal substance derived from human blood plasma) and gentamicin (10 mg/L; medicinal substance).	
QUALITY CONTROL	
<ul style="list-style-type: none">pH between 7.30-7.90 (Release criteria: 7.30-7.60) Osmolality: 270-290 mOsm/kg Endotoxins (USP <85>): < 0.25EU/ml Sterility test by the current Ph. Eur. 2.6.1/ USP <71>: > No growth Human sperm survival assay (% motility compared with control after 24 hours): ≥ 80% One-cell mouse embryo assay (% blastocysts after 96 hours, exposure time to test medium: 60 minutes): ≥ 80% Chemical composition Use of Ph Eur or USP grade products if applicable Certificate of analysis and MSDS are available upon request	
SpermWash® is sterilized by aseptic processing techniques.	
PRECAUTIONS AND WARNINGS	
<ul style="list-style-type: none">Aseptic technique should be used to avoid possible contamination, even when the products contain gentamicin. Always wear protective clothing when handling specimens. All blood products should be treated as potentially infectious. Source material used to manufacture this product was tested and found non-reactive for HbsAg and negative for Anti-HIV-1/2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents. Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include: select and monitor donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/ removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis. Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to Gynotec B.V. and, if applicable, to the competent authority of the EU Member State in which the user and/or patient is established.	
PRE-USE CHECKS	
<ul style="list-style-type: none">Do not use the product if the seal of the container is opened or defect when the product is delivered. Do no use if the products shows any evidence of microbial contamination or becomes cloudy. Do not use after expiry date. Do not freeze before use. Do not re-sterilize after opening. Keep in its original packaging until the day of use. Depending on the number of procedures that will be performed on one day, remove the required volume of medium under aseptic conditions in an appropriate sterile recipient. This is in order to avoid multiple openings/warming cycles of the medium. Discard excess (unused) media. Gentamicin must not be used on a patient that has a known allergy to gentamicin or similar antibiotics.	
STORAGE CONDITIONS	
<ul style="list-style-type: none">Store products between 2-8 °C. Keep away from (sun)light. After opening the container, do not use the product longer than 7 days. Sterile conditions must be maintained and product must be stored at 2°C – 8°C Discard the devices in accordance with local regulations for disposal of medical devices.	
INSTRUCTIONS FOR USE	
Method	
As this medium does contain HEPES, the medium should not be placed in a CO ₂ incubator prior to use. Pre-incubation in incubator (without CO ₂) or heating plate is required to ensure that the temperature of the medium is 37 °C at use. Make sure that the lid is closed when bottles are placed in an incubator with CO ₂ as this will lower the pH to 6.90-6.95.	
Each laboratory should consult its own validated procedures, optimized for its individual medical program. When performing IUI or ET, follow the instructions of the specific catheters used. In addition, some additional information on the following procedures is provided below:	
Washing of spermatozoa (suggested procedure)	
The washing of spermatozoa can be done at room temperature or at 37 °C.	
<ol style="list-style-type: none">Add 5ml SpermWash® to the native semen sample and mix. Centrifuge for 5-10 minutes at approximately 300g. Remove supernatant and leave about 0.5ml of semen in the centrifuge tube. Add 5ml SpermWash® to the test tube. Mix the solution gently until the pellet is completely dissolved. Centrifuge again for 5-10 minutes at 300g. Resuspend in a suitable volume of SpermWash®.	
Swim-up procedure (according to WHO, 2021)	
<ol style="list-style-type: none">Gently layer 1.5ml SpermWash® over 1ml of semen in a conical based centrifuge tube. Incline the tube at an angle of 45° and incubate for 1 hour at 37 °C. Gently turn to the upright position and remove the uppermost 1ml. Dilute this aliquot of motile cells with 2-5ml SpermWash®. Centrifuge for 5-10 minutes at 300g, remove supernatant and finally resuspend in 0.5ml of SpermWash®.	
SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)	
The SSCP for SpermWash® describes safety and performance characteristics for the media and is available on the website of Gynotec B.V. (www.gynotec.nl).	
For further questions regarding to the safety and performance, please contact Gynotec B.V. for customer or technical support.	

Catalogue number	<i>SW-0002</i>	<i>25 x 2 ml SpermWash®</i>	Customer-technical support
<i>Catalogusnummer</i>	<i>SW-0010</i>	<i>25 x 10 ml SpermWash®</i>	<i>Klanten-technische ondersteuning</i>
<i>Référence catalogue</i>	<i>SW-0050</i>	<i>1 x 50 ml SpermWash®</i>	<i>Support clients-support technique</i>
<i>Bestellnummer</i>	<i>SW-0100</i>	<i>1 x 100 ml SpermWash®</i>	<i>Kundendienst-technischer support</i>
<i>Número de catálogo</i>	<i>SW-0250</i>	<i>1 x 250 ml SpermWash®</i>	<i>Atención al cliente-asistencia técnica</i>
<i>Numero di catalogo</i>	<i>SW-0500</i>	<i>1 x 500 ml SpermWash®</i>	<i>Assistenza clienti-supporto tecnico</i>
<i>Número de catálogo</i>			<i>Αποίο αλ cliente-técnico</i>
<i>Αριθμός καταλόγου</i>			<i>Υποστηριξη πελατων-τεχνική υποστηριξη</i>
<i>Katalog numaralari</i>			<i>Müsteri desteği/teknik destek</i>

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