SpermFilter® Stock Solution 100%

ID: FP09 I13_SFSS100 R01 A.1
Date: 01/08/2023

Catalogue number

| SpermFilter® Stock Solution 100% | |
|----------------------------------|---|
| SPF100-0010 | 5 x 10 ml SpermFilter® Stock Solution 100% |
| SPF100-0050 | 1 x 50 ml SpermFilter® Stock Solution 100% |
| SPF100-0100 | 1 x 100 ml SpermFilter® Stock Solution 100% |
| SPF100-0250 | 1 x 250 ml SpermFilter® Stock Solution 100% |
| SPF100-0500 | 1 x 500 ml SpermFilter® Stock Solution 100% |

APPLICATION

SpermFilter® Stock Solution 100% is a stock solution for semen preparation. It is an isotonic gradient for semen preparation with a density of approximately 1.12 g/ml. SpermFilter® Stock Solution 100% can be used for semen preparation for Intra Uterine Insemination (IUI), In Vitro Fertilization (IVF) and IntraCytoplasmic Sperm Injection (ICSI).

For professional use only.

COMPOSITION

SpermFilter® Stock Solution 100% consists of silane-coated colloidal silica particles suspended in HEPES-buffered EBSS (Earle's balanced salt solution).

QUALITY CONTROL

- pH between 7.20-7.90 (Release criteria: 7.20-7.60)
- Osmolality: 300-330 mOsm/kg
- Density: 1.1150-1.1250 g/ml
- Endotoxins (USP <85>): < 0.5EU/ml
- Sterility test according to the current Ph. Eur. 2.6.1/USP <71>:
No growth
- Human sperm survival assay:
 - ◇ % motility compared with control after 4 hours: ≥ 80%
 - ◇ % motility compared with control after 24 hours: ≥ 75%
- Chemical composition
- Use of Ph Eur or USP grade products if applicable
- Not MEA tested
- Certificate of analysis and MSDS are available upon request

SpermFilter® Stock Solution 100% is sterilized by aseptic processing techniques.

PRECAUTIONS AND WARNINGS

- Aseptic technique should be used to avoid possible contamination.
- Always wear protective clothing when handling specimens.
- All human, organic material should be considered potentially infectious. Handle all specimens as if capable of transmitting HIV or hepatitis.
- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to Gynotec B.V. and, if applicable, to the competent authority of the EU Member State in which the user and/or patient is established.

PRE-USE CHECKS

- Do not use the product if the seal of the container is opened or defect when the product is delivered.
- Do not use if the products shows any evidence of microbial contamination or becomes cloudy.
- Do not use after expiry date.
- Do not freeze before use.
- Do not re-sterilize after opening.
- Keep in its original packaging until the day of use.
- Depending on the number of procedures that will be performed on one day, remove the required volume of medium under aseptic conditions in an appropriate sterile recipient. This is in order to avoid multiple openings/warming cycles of the medium. Discard excess (unused) media.

STORAGE CONDITIONS

- Store products between 2-25 °C. Once opened: store between 2-8°C.
- Keep away from (sun)light.
- After opening the container, do not use the product longer than 7 days. Sterile conditions must be maintained and product must be stored at 2-8°C.
- Discard the devices in accordance with local regulations for disposal of medical devices.

CALCULATIONS OF G-FORCES

The g-force of your centrifuge can be calculated using this formula:

$$g = 1.118 \times r \times \text{rpm}^2 \text{ OR } \text{rpm} = \text{Square root } \{g / 1.118 \times r\}$$

r = radius of centrifuge in mm

rpm = rotation per minute / 1000

Example 1

r = 100 mm

rpm = 1800 rotations per minute

$$g = 1.118 \times 100 \times 3.24 = 362 \text{ g}$$

Example 2

r = 100 mm g = 350g

$$\text{rpm} = \text{SQR} \{350 / (1.118 \times 100)\} = 1.77$$

rpm = 1770 rotations per minute

INSTRUCTIONS FOR USE

Method

Each laboratory should establish its own validated procedures and protocols.

Instructions for preparation of gradients

Mix the density gradient bottles by 5 bottle inversions before use.

We advise to prepare a dual gradient system (45%– 90% or 40% – 80%) starting from SpermFilter® Stock Solution 100%. If preferred a multi-layer can be used as well (e.g. 45% – 70% – 90%). To prepare a 90% gradient, mix 1 part SpermTec® Wash or SpermWash® (distributed by Gynotec B.V.) with 9 parts SpermFilter® Stock Solution 100%. A 45% gradient is prepared by mixing 5.5 parts SpermTec® Wash or SpermWash® (distributed by Gynotec B.V.) with 4.5 parts SpermFilter® Stock Solution 100%. Alternatively any HEPES-buffered EBSS-based medium can be used for the preparation of the gradients.

Note: Gradients should be prepared and repacked under sterile conditions (e.g. LAF bench ISO Class 5). For optimal results, prepare the gradients maximum 24 hours prior to use, store at 2-8 °C and warm gradients to room temperature or 37 °C one hour before use. Mix well after diluting the SpermFilter® Stock Solution 100% solution.

Instructions for sperm selection with fresh semen samples: example using a 45%-90% gradient system (but other gradients are possible)

- Transfer 2.5ml of the prepared 45% gradient into a sterile disposable centrifuge tube.
- Place 2.5ml of the prepared 90% gradient under the prepared 45% gradient. Take care that the two layers are distinctly separated. This is done by placing the tip of the pipette on the bottom of the test tube and slowly dispensing the 90% gradient. This two layer gradient is stable for up to two hours.
- Gently place up to 2.5ml of liquefied semen onto the 45% gradient using a transfer pipette. Do not use a higher volume than the volume of the individual gradient layers or more than 10⁹ cells.
- Centrifuge for 15 to 18 minutes at 350g to 400g. When this centrifugation is completed you may not be able to visibly see a pellet. If so, it is essential to continue the procedure with a second centrifugation of 3 to 5 minutes.
- Remove supernatant down to the pellet.
- Add 2-3ml of sperm washing medium (e.g. SpermTec® Wash or SpermWash® distributed by Gynotec B.V.) and resuspend the pellet.
- Centrifuge for 8 to 10 minutes at 300g. Higher sperm concentration will require the maximum 10 minutes centrifugation to ensure a complete and thorough sperm wash.
- Remove supernatant down to the pellet and repeat steps 6 and 7.
- Remove supernatant and replace with a suitable volume of appropriate medium.

Instructions for sperm selection with frozen semen samples: example using a 45%-90% gradient system (but other gradients are possible)

1. Transfer 1ml of the prepared 45% gradient into a sterile disposable centrifuge tube.
2. Place 1ml of the prepared 90% gradient under the prepared 45% gradient. Take care that the two layers are distinctly separated. This is done by placing the tip of the pipette on the bottom of the test tube and slowly dispensing the 90% gradient. This two-layer gradient is stable for up to two hours.
3. Gently place the thawed semen sample onto the 45% gradient using a transfer pipette (0.5-1ml).
4. Centrifuge for 15-20 minutes at 350g.
5. Remove supernatant down to no less than the 0.5ml mark above the pellet.
6. Add 2-3ml of sperm washing medium (e.g. SpermTec® Wash or SpermWash® distributed by Gynotec B.V.) and resuspend the pellet.
7. Centrifuge for 8 to 10 minutes at 300g.
8. Remove supernatant down to the pellet and repeat steps 6 and 7.
9. Remove supernatant and replace with a suitable volume of appropriate medium.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

The SSCP for SpermFilter® Stock Solution 100% describes safety and performance characteristics for the media and is available on the website of Gynotec B.V. (www.gynotec.nl).

For further questions regarding to the safety and performance, please contact Gynotec B.V. for customer or technical support.

SYMBOL GLOSSARY

| Symbols as defined in ISO 15223-1 | | | |
|--|---|---|---|
|  | Catalogue number |  | Medical device |
|  | Use-by date |  | Cautions |
|  | Keep away from sunlight |  | Sterilized using aseptic processing techniques |
|  | Consult instructions for use |  | Do not use if package is damaged and consult instructions for use |
|  | Date of manufacture |  | Single sterile barrier system |
|  | Temperature limit |  | Single sterile barrier system with protective packaging outside |
|  | Batch code |  | Distributor |
|  | Manufacturer |  | Do not resterilize |
| Symbol as defined in EU Commission MDR 2017/745 | | | |
|  | CE marking by Notified Body 2797 | | |
| Symbol as defined in the Resolution of the Cabinet of Ministers of Ukraine No. 1184 dated 30.12.2015 | | | |
|  | Mark of conformity with technical regulations accompanied by the identification number of the designated conformity assessment body | | |

CUSTOMER-TECHNICAL SUPPORT

 **Gynotec B.V.**
 Jonckherenhof 7
 6581 GC Malden The Netherlands
 Tel: +31 (0)243 58 65 82
 Fax +31 (0)243 58 13 55
 URL: <http://www.gynotec.nl>
 E-mail: info@gynotec.nl

 **FertiPro NV**
 Industriepark Noord 32
 8730 Beernem Belgium
 Tel: +32 (0)50 79 18 05
 Fax: +32 (0)50 79 17 99
 URL: <https://www.fertipro.com>
 E-mail: info@fertipro.com

UKRAINIAN AUTHORIZED REPRESENTATIVE

«Stoik Company» LLC
 C 04136 ap. A-612, building A, 13 Ivana Vygovskogo str., Kyiv, Ukraine
 Tel: +38 (0)67 112-92-92
 URL: www.stoik.com.ua
 E-mail: stoik.ivf@gmail.com

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