♦ APPLICATION

SpermTec® Cryo is a four times concentrated ready-to-use HEPES buffered cryopreservation medium for freezing human sperm for further use in Assisted Reproductive Technologies (ART).

For professional use only.

COMPOSITION
SpermTec® Cryo is a HEPES buffered cryopreservation medium which contains physiologic salts, glycine, glucose, lactate and the cryoprotectants glycerol (27%), sucrose and Human Serum Albumin (4.0 g/Liter; medicinal substance derived from human blood plasma) in order to protect the sperm from damage during the freezing procedure.

MATERIAL NOT INCLUDED

- Sperm freezing straws (e.g. CBS™ high security sperm straws) or cryo tubes
- Freezing tank with liquid nitrogen

QUALITY CONTROL

- pH between 7.20 – 7.90 (Release criteria: 7.20 – 7.60)
- Endotoxins (USP <85>): < 0.25EU/ml
- Sterility test by the current Ph. Eur. 2.6.1 / USP <71>: No growth
- Human sperm survival assay (% motility compared with control after 4 hours): ≥ 80%
- Chemical composition
- Use of Ph Eur or USP grade products if applicable
- Notation and testing
- Certificate of analysis and MSDS are available upon request

SpermTec® Cryo is sterilized by aseptic processing techniques.

PRECAUTIONS AND WARNINGS

- Aseptic techniques should be used to avoid possible contamination. SpermTec® Cryo does not contain any antibiotics.
- Always wear protective clothing when handling specimens.
- All blood products should be treated as potentially infectious. Source reduction studies have shown this product was tested and found non-reactive for HbsAg and negative for Anti-HIV-1/-2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents.
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps (e.g. Sterilization/Freezing).
- Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, human blood specimens as if capable of transmitting HIV or hepatitis.
- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to Gynotec B.V. and, if applicable, the competent authority of the EU Member State in which the user and/or patient is established.

PRE-USE CHECKS

- Do not use the product if the seal of the container is opened or defect when the product is delivered.
- Do not use if the products shows any evidence of microbial contamination or becomes cloudy.
- Do not use if the expiry date has expired.
- Do not freeze before use.
- Do not re-sterilize after opening.
- Keep in its original packaging until the day of use.
- Depending on the number of procedures that will be performed on one day, remove the required volume of medium under aseptic conditions in an appropriate sterile recipient. This is in order to avoid multiple openings/warming cycles of the medium. Discard excess (unused) media.

STORAGE CONDITIONS

- Store refrigerated 2-8 °C.
- Keep away from (sun)light.
- After opening the container, do not use the product longer than 7 days. Sterile conditions are maintained and product must be stored at 2-8 °C.
- Discard the devices in accordance with local regulations for disposal of medical devices.

INSTRUCTIONS FOR USE

Sperm preparation

Before freezing
Cryopreservation is ideally performed on native semen samples.

In case of very low sperm concentrations it is advisable to concentrate the sperm before freezing. In case of very high sperm concentration, sperm can be diluted with a sperm washing medium before freezing (e.g. SpermWash® / SpermTec® Wash distributed by Gynotec B.V.), if preferred.

After freezing

In case needed, use sperm preparation techniques after thawing the semen to eliminate dead sperm cells and debris. Dilute the concentrated sperm in a suitable insemination medium.

Method

Ensure that all media are well mixed and at room temperature before use.

Freezing

- Allow the semen to liquefy at room temperature for 30 minutes.
- Add 1 part SpermTec® Cryo to 3 parts semen: add the medium in drops while gently swirling.
- Cautions: to avoid cold-shock, make sure SpermTec® Cryo is at room temperature.
- Leave the mixture for 10 minutes at room temperature for 10 minutes.
- Suck the sample/medium mixture into the freezing straws, leaving approximately 1.5cm of air at the end of the straw.
- Make sure that the outer surface of the straw is clean (no sperm) and seal the straw.
- Dry, if applicable each straw with a dust free cloth and shake to move the air-bubble to the center of the straw.
- Freeze the straws in vapor phase by placing the straws horizontally (e.g. on a Styrofoam board) in a liquid nitrogen bath for (at least) 15 minutes.
- Quickly transfer the straws into liquid nitrogen and store at -196 °C.

Thawing

- Retrieve the required straws from the liquid nitrogen.
- Place the straws in tap water for 5 minutes (at room temperature or 37 °C).
- Cut off the end of the straw, place the open end inside a container (e.g. test tube) and tap the straw against the side of the container to allow complete evacuation of the mixture.
- Dilute the concentrated sperm in a suitable insemination medium (e.g. SpermWash® / SpermTec® Wash) (at least 3ml per 0.5ml semen) and mix thoroughly.
- Centrifuge at 300-350g for 15 minutes.
- Resuspend pellet in a suitable insemination medium (e.g. SpermWash® / SpermTec® Wash) and access recovery.

Notes: For additional details on the method, each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for its individual medical program.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)
Although damaging effects on sperm motility, vitality, DNA fragmentation and fertilizing potential upon thawing are inherent to the cryopreservation process, these remain within acceptable limits for SpermTec® Cryo as reported in scientific literature. The SSCP for SpermTec® Cryo describes safety and performance characteristics for the media and is available on the website of Gynotec B.V. (www.gynotec.nl).

For further questions regarding to the safety and performance, please contact Gynotec B.V. for customer or technical support.

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