INTENDED USE

SpermFilter® Ready-to-use gradients (80% and 45%) are used as sperm preparation method for intrauterine insemination (IUI), in vitro fertilization (IVF), intracytoplasmic sperm injection (ICSI), and related assisted reproductive technologies (ART).

INDICATIONS FOR USE

For use during ART procedures of patients and couples undergoing infertility treatments.

INTENDED USER

The intended users are ART professionals (lab technicians, embryologists, or medical doctors).

TARGET PATIENT POPULATION

The target patient population consists of patients and couples undergoing infertility treatments.

COMPOSITION

SpermFilter® Ready-to-use gradients consist of silane-coated colloidal silica particles suspended in HEPES-buffered EBSS (Earle's balanced salt solution) supplemented with phenol red, 0.4-2.2 g/l Human Serum Albumin (medicinal substance derived from human blood plasma) and 10mg/l gentamicin (medicinal substance).

QUALITY CONTROL

- pH between 7.20-7.90 (Release criteria: 7.20-7.60)
- Osmolality:
 - 80% gradient: 290-330 mOsm/kg
 - 45% gradient: 280-310 mOsm/kg
- Density:
 - 80% gradient: 1.097-1.107 g/ml
 - 45% gradient: 1.0500-1.0700 g/ml
- Sterility test according to the current Ph. Eur. 2.6.1/ USP <71>: No growth
- Human sperm survival assay:
 - % motility compared with control after 4 hours: > 75%
 - % motility compared with control after 24 hours: > 75%
- Chemical composition
 - Use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/elimination of viruses, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of previous viruses transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.
 - Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FertiPro NV and, if applicable, to the competent authorities of the EU member state in which the user and/or patient is established.

Sterilized by using aseptic processing techniques.

PRECAUTIONS AND WARNINGS

- Aseptic technique should be used to avoid possible contamination, even when the products contain gentamicin.
- Always wear protective clothing when handling specimens.
- All blood products should be treated as potentially infectious. Source material used in the manufacture of this product was tested and found non-reactive for HbsAg and negative for Anti-HIV-1/-2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non-reactive. No known test method can offer assurances that products derived from human blood will not transmit infectious agents.
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/elimination of viruses, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of previous viruses transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.
- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FertiPro NV and, if applicable, to the competent authorities of the EU member state in which the user and/or patient is established.

PRE-USE CHECKS

- Do not use the product if the seal of the container is broken or if the product is damaged.
- Do not use if the products shows any evidence of microbial contamination or becomes cloudy.
- Do not use after expiry date.
- Do not freeze before use.
- Do not re-sterilize after opening.
- Keep in its original packaging until the day of use.
- Depending on the number of procedures that will be performed on one day, remove the required volume of medium under aseptic conditions in an appropriate sterile recipient. This is in order to avoid multiple openings/warming cycles of the medium. Discard excess (unused) media.
- Gentamicin should not be used on a patient that has a known allergy to gentamicin or similar antibiotics.

STORAGE CONDITIONS

- Store products between 2-8 °C.
- Keep away from direct sunlight.
- After opening the container, do not use the product longer than 7 days. Sterile conditions must be maintained and product must be stored at 2 °C – 8 °C.
- Discard the devices in accordance with local regulations for disposal of medical devices.

CALCULATIONS OF G-FORCES

The g-force of your centrifuge can be calculated using this formula:

$$g = \frac{r \times \text{rpm}^2}{1.118 \times r \times \text{rpm}^2} \quad \text{OR} \quad \text{rpm} = \sqrt{\frac{\text{Square root}}{(g / 1.118 \times r)}}$$

Example 1 r = 100 mm rpm = 1800 rotations per minute g = 1.118 x 100 x 3.24 = 362 g	Example 2 r = 100 mm g = 350g rpm = $\sqrt{350 / (1.118 \times 100)}$ = 1.77 rpm = 1770 rotations per minute
---	---

INSTRUCTIONS FOR USE

Method

Each laboratory should establish its own validated procedures and protocols.

Fresh semen samples

Mix the density gradient bottles by 5 bottle inversions before use.

- Bring all components of the system and samples to room temperature or to 37 °C.
- Transfer 2.5ml of SpermFilter® Ready-to-use gradient 45% into a sterile disposable centrifuge tube.
- Place 2.5ml of SpermFilter® Ready-to-use gradient 80% under the SpermFilter® Ready-to-use gradient 45% layer. Take care that the two layers are distinctly separated. This is done by placing the tip of the pipette on the bottom of the test tube and slowly dispensing the SpermFilter® Ready-to-use gradient 80%. This two layer gradient is stable for up to two hours.
- Gently place up to 2.5ml of liquefied semen onto the SpermFilter® Ready-to-use gradient 45% using a transfer pipette. Do not use a higher volume than the volume of the individual gradient layers or more than 4x10⁶ cells.
- Centrifuge for 15 to 18 minutes at 350g to 400g. When this centrifugation is completed you may not be able to visibly see a pellet. If so, it is essential to continue the procedure with a second centrifugation of 3 to 5 minutes.
- Remove supernatant down to the pellet.
- Add 2-3ml of sperm washing medium (e.g. SpermTec® Wash or SpermWash® distributed by Gynotec B.V.) and resuspend the pellet.
- Centrifuge for 8 to 10 minutes at 300g. Higher sperm concentration will require the maximum 10 minutes centrifugation to ensure a complete and thorough sperm wash.
- Remove supernatant down to the pellet and repeat steps 7 and 8.
- Remove supernatant and replace with a suitable volume of appropriate media.

Frozen semen samples

Mix the density gradient bottles by 5 bottle inversions before use.

- Bring all components of the system and samples to room temperature or to 37 °C.
- Transfer 1ml of SpermFilter® Ready-to-use gradient 45% into a sterile disposable centrifuge tube.
- Place 1ml of SpermFilter® Ready-to-use gradient 80% under the SpermFilter® Ready-to-use gradient 45%. Make sure that the two layers are distinctly separated. This is done by placing the tip of the pipette on the bottom of the test tube and slowly dispensing the SpermFilter® Ready-to-use gradient 80%. This two layer gradient is stable for up to two hours.
- Gently place the thawed semen sample onto the SpermFilter® Ready-to-use gradient 45% using a transfer pipette (0.5ml – 1ml).
- Centrifuge for 15-20 minutes at 350g.
- Remove supernatant down to no less than the 0.5ml mark above the pellet.
- Add 2-3ml of sperm washing medium (e.g. SpermTec® Wash or SpermWash® distributed by Gynotec B.V.) and resuspend the pellet.
- Centrifuge for 8 to 10 minutes at 300g.
- Remove supernatant down to the pellet and repeat steps 7 and 8.
- Remove supernatant and replace with a suitable volume of appropriate media.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

The SSCP for SpermFilter® Ready-to-use gradient 80% and 45% describes safety and performance characteristics for the media and is available on the website of Gynotec B.V. (www.gynotec.nl).

For further questions regarding to the safety and performance, please contact Gynotec B.V. for customer or technical support.

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- Bring all components of the system and samples to room temperature or to 37 °C.
- Transfer 1ml SpermFilter® Ready-to-use gradient 45% into a sterile well-washed centrifuge tube.
- Place 1ml SpermFilter® Ready-to-use gradient 80% under the SpermFilter® Ready-to-use gradient 45%. Make sure that the two layers are distinctly separated. This is done by placing the tip of the pipette on the bottom of the test tube and slowly dispensing the SpermFilter® Ready-to-use gradient 80%. This two layer gradient is stable for up to two hours.
- Gently place the thawed sperm sample onto the SpermFilter® Ready-to-use gradient 45% using a transfer pipette (0.5 - 1ml).
- Centrifuge for 15-20 minutes at 350g.
- Remove supernatant down to no less than the 0.5ml mark above the pellet.
- Add 2-3ml of sperm washing medium (e.g. SpermTec® Wash or SpermWash® distributed by Gynotec B.V.) and resuspend the pellet.
- Centrifuge for 8 to 10 minutes at 300g.
- Remove supernatant down to the pellet and repeat steps 7 and 8.
- Remove supernatant and replace with a suitable volume of appropriate media.

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- Bring all components of the system and samples to room temperature or to 37 °C.
- Transfer 1ml SpermFilter® Ready-to-use gradient 45% into a sterile disposable centrifuge tube.
- Place 1ml SpermFilter® Ready-to-use gradient 80% under the SpermFilter® Ready-to-use gradient 45%. Make sure that the two layers are distinctly separated. This is done by placing the tip of the pipette on the bottom of the test tube and slowly dispensing the SpermFilter® Ready-to-use gradient 80%. This two layer gradient is stable for up to two hours.
- Gently place the thawed sperm sample onto the SpermFilter® Ready-to-use gradient 45% using a transfer pipette (0.5 - 1ml).
- Centrifuge for 15-20 minutes at 350g.
- Remove supernatant down to no less than the 0.5ml mark above the pellet.
- Add 2-3ml of sperm washing medium (e.g. SpermTec® Wash or SpermWash® distributed by Gynotec B.V.) and resuspend the pellet.
- Centrifuge for 8 to 10 minutes at 300g.
- Remove supernatant down to the pellet and repeat steps 7 and 8.
- Remove supernatant and replace with a suitable volume of appropriate media.

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