IVF Basics® HTF IVF Basics® HTF 0.4% HSA

IVF Basics® HTF HEPES IVF Basics® HTF HEPES 0.4% HSA

EN

USED ABBREVIATIONS

ICSI Intracytoplasmatic Sperm Injection
IVF In Vitro Fertilization
IUI Intra uterine insemination
HSA Human Serum Albumin
LAL Limulus Amebocyte Lysate
SSCP Summary of Safety and Clinical Performance
MSDS Material Safety Data Sheet

GENERAL

IVF Basics® HTF (0.4% HSA) and IVF Basics® HTF HEPES (0.4% HSA) have been developed for washing, handling, conditioning and culturing of human gametes and embryos.
► IVF Basics® HTF (0.4% HSA) should be pre-equilibrated overnight or at least for a minimum of 4 hours before use at 37°C and 5% CO₂.
► IVF Basics® HTF HEPES (0.4% HSA) contains HEPES, no CO₂ incubation is required.

APPLICATION

Intended use
IVF Basics® HTF / IVF Basics® HTF 0.4% HSA are used for washing and handling of human gametes, and for handling and conditioning of human embryos. The media are used for the following procedures:

- Washing of human ova
- Washing/handling of spermatozoa
- Swim-up of spermatozoa
- Production of density gradient media
- Washing/handling for human embryos
- IUI, IVF, ICSI
- Embryo culture from day 1 to blastocyst stage
- Embryo transfer

IVF Basics® HTF HEPES / IVF Basics® HTF HEPES 0.4% HSA are used for washing and handling of human gametes, and for handling and conditioning of human embryos. The media are used for the following procedures:

- Flushing for oocyte collection
- Washing/handling of human ova
- Washing/handling of spermatozoa
- Swim-up of spermatozoa
- Production of density gradient media
- Washing/handling for human embryos
- IUI, ICSI
- Embryo transfer

Indications for use

For use during ART procedures of patients and couples undergoing infertility treatments.

Intended user

The intended users are ART professionals (lab technicians, embryologists, or medical doctors).

Target patient population

The target patient population consists of patients and couples undergoing infertility treatments.

COMPOSITION

IVF Basics® HTF and IVF Basics® HTF HEPES consist of a balanced physiologic salt solution, supplemented with carbohydrate energy sources such as glucose, pyruvate and lactate, and 10mg/liter gentamicin sulphate (medical substance). IVF Basics® HTF is bicarbonate buffered, while IVF Basics® HTF HEPES additionally contains HEPES.

It is strongly suggested to supplement IVF Basics® HTF and IVF Basics® HTF HEPES with 4g/liter HSA (medical substance derived from human blood plasma) before use.

IVF Basics® HTF 0.4% HSA and IVF Basics® HTF HEPES 0.4% HSA are pre-equilibrated at 37°C and 5% CO₂. IVF Basics® HTF HEPES 0.4% HSA is ready-to-use.

The MSDS for IVF Basics® HTF media also describes the product composition and is available on the website of Gynotec B.V.

MATERIAL NOT INCLUDED

► IVF Basics® HTF and IVF Basics® HTF HEPES (media without HSA): pharmaceutical grade Human Albumin Solution with registered Plasma Master File

QUALITY CONTROL

- pH: 7.20-7.60 (At 5% CO₂ for IVF Basics® HTF (with HSA))
- Osmolality: 270-290 mOsm/kg
- Osmolality: 270-290 mOsm/kg
- Sterility test by the current Ph. Eur. 2.6.1/USP <71>: No growth
- One-cell mouse embryo assay (blastocysts after 96h): ≥80%
- Human sperm survival assay (5% motility compared with control after 24 hours): ≥80%
- Chemical composition
- Use of Ph. Eur. or USP grade products if applicable
- Certificate of analysis and MSDS are available upon request

Sterilized by using aseptic processing techniques.

PRECAUTIONS AND WARNINGS

- Aseptic techniques should be used to avoid possible contamination, even when the products contain gentamicin.
- Always wear protective clothing when working with specimens.
- Handle specimens as if capable of transmitting HIV or hepatitis.
- All blood products should be treated as potentially infectious. Source material used to manufacture this product was tested and found non-reactive for HbsAg and negative for Anti-HIV-1/2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non-evaluated. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents.
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include detection of donors, screening of individual donations and plasma

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► Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FeriPro NV and, if applicable, to the competent authority of the EU Member State in which the user and/or patient is established.

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IVF Basics® HTF

Product code	Product description
HTF-B-10	1 x 10 ml of IVF Basics® HTF
HTF-B-30	1 x 30 ml of IVF Basics® HTF
HTF-B-60	1 x 60 ml of IVF Basics® HTF
HTF-B-120	1 x 120 ml of IVF Basics® HTF
HTF-B-250	1 x 250 ml of IVF Basics® HTF
HTF-B-500	1 x 500 ml of IVF Basics® HTF

IVF Basics® HTF HEPES

Product code	Product description
HTF-B-10	1 x 10 ml of IVF Basics® HTF
HTF-B-30	1 x 30 ml of IVF Basics® HTF
HTF-B-60	1 x 60 ml of IVF Basics® HTF
HTF-B-120	1 x 120 ml of IVF Basics® HTF
HTF-B-250	1 x 250 ml of IVF Basics® HTF
HTF-B-500	1 x 500 ml of IVF Basics® HTF

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