SpermWash®

ID: FP09 I08_SW R01 A.3
Date: 16/03/2026

Catalogue number

SpermWash®	
SW-0002	25 x 2 ml SpermWash®
SW-0010	25 x 10 ml SpermWash®
SW-0050	1 x 50 ml SpermWash®
SW-0100	1 x 100 ml SpermWash®
SW-0250	1 x 250 ml SpermWash®
SW-0500	1 x 500 ml SpermWash®

APPLICATION

Intended use

SpermWash® is used for washing of human gametes (sperm and oocytes), sperm swim-up procedures, intrauterine insemination (IUI) of the spermatozoa and intracytoplasmic sperm injection (ICSI). SpermWash® is also used for human embryo washing and holding, and for embryo transfer in the uterine cavity.

Indications for use

For use during ART procedures of patients and couples undergoing infertility treatments.

Intended user

The intended users are ART professionals (lab technicians, embryologists, or medical doctors).

Target patient population

The target patient population consists of patients and couples undergoing infertility treatments.

COMPOSITION

SpermWash® is a ready-to-use HEPES buffered medium which also contains phenol red, bicarbonate, physiologic salts, glucose, lactate, pyruvate and Human Serum Albumin (4.0 g/l; medicinal substance derived from human blood plasma) and gentamicin (10 mg/l, medicinal substance).

QUALITY CONTROL

- pH between 7.30-7.90 (Release criteria: 7.30-7.60)
- Osmolality: 270-290 mOsm/kg
- Endotoxins (USP <85>): < 0.25EU/ml
- Sterility test by the current Ph. Eur. 2.6.1/ USP <71>: No growth
- Human sperm survival assay (% motility compared with control after 24 hours): ≥ 80%
- One-cell mouse embryo assay (% blastocysts after 96 hours, exposure time to test medium: 60 minutes): ≥ 80%
- Chemical composition
- Use of Ph Eur or USP grade products if applicable
- Certificate of analysis and MSDS are available upon request

SpermWash® is sterilized by aseptic processing techniques.

PRECAUTIONS AND WARNINGS

- Aseptic technique should be used to avoid possible contamination, even when the products contain gentamicin.
- Always wear protective clothing when handling specimens.
- All blood products should be treated as potentially infectious. Source material used to manufacture this product was tested and found non-reactive for HbsAg and negative for Anti-HIV-1/-2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be nonelevated. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents.
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses

and other pathogens. There are no reports of proven viruses transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FertiPro NV and, if applicable, to the competent authority of the EU Member State in which the user and/or patient is established.

PRE-USE CHECKS

- Do not use the product if the seal of the container is opened or defect when the product is delivered.
- Do not use if the products shows any evidence of microbial contamination or becomes cloudy.
- Do not use after expiry date.
- Do not freeze before use.
- Do not re-sterilize after opening.
- Keep in its original packaging until the day of use.
- Depending on the number of procedures that will be performed on one day, remove the required volume of medium under aseptic conditions in an appropriate sterile recipient. This is in order to avoid multiple openings/warming cycles of the medium. Discard excess (unused) media.
- Gentamicin should not be used on a patient that has a known allergy to gentamicin or similar antibiotics.

STORAGE CONDITIONS

- Store products between 2-8 °C.
- Keep away from (sun)light.
- After opening the container, do not use the product longer than 7 days. Sterile conditions must be maintained and product must be stored at 2°C – 8°C.
- Discard the devices in accordance with local regulations for disposal of medical devices.

INSTRUCTIONS FOR USE

Method

As this medium does contain HEPES, the medium should not be placed in a CO₂ incubator prior to use. Preincubation in incubator (without CO₂) or heating plate is required to ensure that the temperature of the medium is 37 °C at use. Make sure that the lid is closed when bottles are placed in an incubator with CO₂ as this will lower the pH to 6.90-6.95.

Each laboratory should consult its own validated procedures, optimized for its individual medical program. When performing IUI or ET, follow the instructions of the specific catheters used. In addition, some additional information on the following procedures is provided below:

Washing of spermatozoa (suggested procedure)

The washing of spermatozoa can be done at room temperature or at 37 °C.

- Add 5ml SpermWash® to the native semen sample and mix. Centrifuge for 5-10 minutes at approximately 300g.
- Remove supernatant and leave about 0.5ml of semen in the centrifuge tube.
- Add 5ml SpermWash® to the test tube. Mix the solution gently until the pellet is completely dissolved.
- Centrifuge again for 5-10 minutes at 300g.
- Resuspend in a suitable volume of SpermWash®.

Swim-up procedure (according to WHO, 2021)

- Gently layer 1.5ml SpermWash® over 1ml of semen in a conical based centrifuge tube.
- Incline the tube at an angle of 45° and incubate for 1 hour at 37 °C.
- Gently turn to the upright position and remove the uppermost 1ml.
- Dilute this aliquot of motile cells with 2-5ml SpermWash®. Centrifuge for 5-10 minutes at 300g, remove supernatant and finally resuspend in 0.5ml of SpermWash®.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

The SSCP for SpermWash® describes safety and performance characteristics for the media and is available on the website of Gynotec B.V. (www.gynotec.nl).

For further questions regarding to the safety and performance, please contact Gynotec B.V. for customer or technical support.

SYMBOL GLOSSARY

Symbols as defined in ISO 15223-1			
	Catalogue number		Cautions
	Use-by date		Sterilized using aseptic processing techniques
	Keep away from sunlight		Do not use if package is damaged and consult instructions for use
	Consult instructions for use		Single sterile barrier system
	Date of manufacture		Single sterile barrier system with protective packaging outside
	Temperature limit		Distributor
	Batch code		Contains human blood or plasma derivatives
	Manufacturer		Contains a medicinal substance
	Medical device		Do not resterilize
Symbol as defined in EU Commission MDR 2017/745			
	CE marking by Notified Body 2797		
Symbol as defined in the Resolution of the Cabinet of Ministers of Ukraine No. 1184 dated 30.12.2015			
	Mark of conformity with technical regulations accompanied by the identification number of the designated conformity assessment body		

CUSTOMER-TECHNICAL SUPPORT

Gynotec B.V.
Rijksweg 150
6581 ET Malden The Netherlands
Tel: +31 (0)243 58 65 82
Fax: +31 (0)243 58 13 55
URL: <http://www.gynotec.nl>
E-mail: info@gynotec.nl

FertiPro NV
Industriepark Noord 32
8730 Beernem Belgium
Tel: +32 (0)50 79 18 05
Fax: +32 (0)50 79 17 99
URL: <https://www.fertipro.com>
E-mail: info@fertipro.com

UKRAINIAN AUTHORIZED REPRESENTATIVE

«Stoik Company» LLC
C 04136 ap. A-612, building A, 13 Ivana Vygovskogo str., Kyiv, Ukraine
Tel: +38 (0)67 112-92-92
URL: www.stoik.com.ua
E-mail: stoik.ivf@gmail.com

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