SpermCryo® All-round

ID: FP09 I11_SPCA R01 A.3
Date: 16/03/2026

Catalogue number

SpermCryo® All-round	
SPCA-0005	25 x 5ml SpermCryo® All-round
SPCA-0010	5 x 10ml SpermCryo® All-round
	25 x 10ml SpermCryo® All-round

APPLICATION

Intended use

SpermCryo® All-round is used for cryopreservation of human sperm for further use in assisted reproductive technologies (ART).

Indications for use

For use during ART procedures of patients and couples undergoing infertility treatments.

Intended user

The intended users are ART professionals (lab technicians, embryologists, or medical doctors).

Target patient population

The target patient population consists of patients and couples undergoing infertility treatments.

COMPOSITION

SpermCryo® All-round is a four times concentrated ready-to-use HEPES buffered cryopreservation medium which contains physiologic salts, glycine, glucose, lactate and the cryoprotectants glycerol (27%), sucrose and Human Serum Albumin (4.0 g/liter; medicinal substance derived from human blood plasma) in order to protect the sperm from damage during the freezing procedure.

MATERIAL NOT INCLUDED

- Sperm freezing straws (e.g. CBS™ high security sperm straws) or cryo tubes
- Freezing tank with liquid nitrogen

QUALITY CONTROL

- pH between 7.20-7.90 (Release criteria: 7.20-7.60)
- Endotoxins (USP <85>): < 0.25EU/ml
- Sterility test by the current Ph. Eur. 2.6.1/ USP <71>: No growth
- Human sperm survival assay (% motility compared with control after 4 hours): ≥ 80%
- Chemical composition
- Use of Ph Eur or USP grade products if applicable
- Not MEA tested
- Certificate of analysis and MSDS are available upon request

SpermCryo® All-round is sterilized by aseptic processing techniques.

PRECAUTIONS AND WARNINGS

- Aseptic technique should be used to avoid possible contamination. SpermCryo® All-round does not contain any antibiotics.
- Always wear protective clothing when handling specimens.
- All blood products should be treated as potentially infectious. Source material used to manufacture this product was tested and found non-reactive for HbsAg and negative for Anti-HIV-1/2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be nonelevated. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents.
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven viruses transmissions with albumin manufactured to European

Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FertiPro NV and, if applicable, to the competent authority of the EU Member State in which the user and/or patient is established.

PRE-USE CHECKS

- Do not use the product if the seal of the container is opened or defect when the product is delivered.
- Do not use if the products shows any evidence of microbial contamination or becomes cloudy.
- Do not use after expiry date.
- Do not freeze before use.
- Do not re-sterilize after opening.
- Keep in its original packaging until the day of use.
- Depending on the number of procedures that will be performed on one day, remove the required volume of medium under aseptic conditions in an appropriate sterile recipient. This is in order to avoid multiple openings/warming cycles of the medium. Discard excess (unused) media.

STORAGE CONDITIONS

- Store refrigerated 2-8 °C.
- Keep away from (sun)light.
- After opening the container, do not use the product longer than 7 days. Sterile conditions must be maintained and product must be stored at 2-8 °C.
- Discard the devices in accordance with local regulations for disposal of medical devices.

INSTRUCTIONS FOR USE

Sperm preparation

Before freezing

Cryopreservation is ideally performed on native semen samples.

In case of very low sperm concentrations it is advisable to concentrate the sperm before freezing. In case of very high sperm concentration, sperm can be diluted with a sperm washing medium before freezing (e.g. SpermWash® / SpermTec® Wash distributed by Gynotec B.V.), if preferred.

After freezing

In case needed, use sperm preparation techniques after thawing the semen to eliminate dead sperm cells and debris. Dilute the concentrated sperm in a suitable insemination medium.

Method

Ensure that all media are well mixed and at room temperature before use.

Freezing

- Allow the semen to liquefy at room temperature for 30 minutes.
- Add 1 part SpermCryo® All-round to 3 parts semen: add the medium in drops while gently swirling.
- Caution: to avoid cold-shock, make sure SpermCryo® All-round is at room temperature.
- Leave the mixture for 10 minutes at room temperature for equilibration.
- Suck the sample/medium mixture into the freezing straws, leaving approximately 1.5cm of air at the end of the straw.
- Make sure that the outer surface of the straw is clean (no sperm) and seal the straw.
- Dry, if applicable each straw with a dust free cloth and shake to move the air-bubble to the center of the straw.
- Freeze the straws in vapor phase by placing the straws horizontally (e.g. on a Styrofoam board) in a liquid nitrogen bath for (at least) 15 minutes.
- Quickly transfer the straws into liquid nitrogen and store at -196 °C.

Thawing

- Retrieve the required straws from the liquid nitrogen.
- Place the straws in tap water for 5 minutes (at room temperature or 37 °C).
- Cut off the end of the straw, place the open end inside a container (e.g. test tube) and tap the straw against the side of the container to allow complete evacuation of the mixture.

- Dilute the concentrated sperm in a suitable insemination medium (e.g. SpermWash® / SpermTec® Wash) (at least 3ml per 0.5ml semen) and mix thoroughly.
- Centrifuge at 300-350g for 15 minutes.
- Resuspend pellet in a suitable insemination medium (e.g. SpermWash® / SpermTec® Wash) and access recovery.

Notes: For additional details on the method, each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for its individual medical program.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

Although damaging effects on sperm motility, vitality, DNA fragmentation and fertilizing potential upon thawing are inherent to the cryopreservation process, these remain within acceptable limits for SpermCryo® All-round as reported in scientific literature. The SSCP for SpermCryo® All-round describes safety and performance characteristics for the media and is available on the website of Gynotec B.V. (www.gynotec.nl).

For further questions regarding to the safety and performance, please contact Gynotec B.V. for customer or technical support.

SYMBOL GLOSSARY

Symbols as defined in ISO 15223-1			
	Catalogue number		Cautions
	Use-by date		Sterilized using aseptic processing techniques
	Keep away from sunlight		Do not use if package is damaged and consult instructions for use
	Consult instructions for use		Single sterile barrier system
	Date of manufacture		Single sterile barrier system with protective packaging outside
	Temperature limit		Distributor
	Batch code		Contains human blood or plasma derivatives
	Manufacturer		Contains a medicinal substance
	Medical device		Do not resterilize
Symbol as defined in EU Commission MDR 2017/745			
	CE marking by Notified Body 2797		
Symbol as defined in the Resolution of the Cabinet of Ministers of Ukraine No. 1184 dated 30.12.2015			
	Mark of conformity with technical regulations accompanied by the identification number of the designated conformity assessment body		

CUSTOMER-TECHNICAL SUPPORT

 **Gynotec B.V.**
 Rijksweg 150
 6581 ET Malden The Netherlands
 Tel: +31 (0)243 58 65 82
 Fax: +31 (0)243 58 13 55
 URL: <http://www.gynotec.nl>
 E-mail: info@gynotec.nl

 **FertiPro NV**
 Industriepark Noord 32
 8730 Beernem Belgium
 Tel: +32 (0)50 79 18 05
 Fax: +32 (0)50 79 17 99
 URL: <https://www.fertipro.com>
 E-mail: info@fertipro.com

UKRAINIAN AUTHORIZED REPRESENTATIVE

«Stoik Company» LLC
 C 04136 ap. A-612, building A, 13 Ivana Vygovskogo str., Kyiv, Ukraine
 Tel: +38 (0)67 112-92-92
 URL: www.stoik.com.ua
 E-mail: stoik.ivf@gmail.com

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